


PCR plates with half skirt – high profile




Product information:

Profile: High profile
 Maximum well volume: 0.2 ml
 Cutaway corner: A12

Features and benefits:

- Highly uniform and thin well walls ensure consistently even and quickest possible transfer of heat. This guarantees reliable and highly reproducible results.
- The raised skirt around each well provides protection against cross contamination and facilitates secure sealing with films and foils, thus protecting against evaporation loss.
- Black alphanumeric labelling facilitates sample identification and traceability during manual filling.
- Production under clean room conditions and independently conducted biological tests allow for outstanding PCR Performance Tested and Biosphere® plus purity certifications.
- Testing of each individual well for 100% leakproofness ensures maximum security for valuable samples.
- Plates can be stacked safely and securely, enabling efficient use of space if storage is limited.



Description	Colour	Purity	Packaging	Order no:
96-well PCR plate with half skirt	Transparent		10 units/bag & 100 units/box	72.1979
96-well PCR plate with half skirt	Transparent		1 unit/blister pack & 20 units/box	72.1979.201
96-well PCR plate with half skirt	White (qPCR-optimised)		10 units/bag & 100 units/box	72.1979.010
96-well PCR plate with half skirt & barcode	Transparent		10 units/bag & 100 units/box	72.1979.003
96-well PCR plate with half skirt, DNA low binding	Transparent		10 units/bag & 100 units/box	72.1979.700
96-well PCR plate with half skirt and flat deck	Transparent		5 units/bag & 100 units/box	72.1979.102
96-well PCR plate with half skirt and flat deck	White (qPCR-optimised)		5 units/bag & 100 units/box	72.1979.132

Other colours and barcoded variants available on request.
 Suitable lid strips and sealing films and foils can be found on pages 24–27.