

# MegaBlock® 96 Well



Free from human DNA,  
DNase/RNase, endotoxins

## MegaBlock® 96 Well

The MegaBlock® 96 Well is ideal for automated processing of sample volumes up to 2.2 ml or handling large quantities of backup samples.

- Wells with alphanumeric coordination system
- Each well 100% leak-proof tested
- Ideal for long-term storage of samples
- Also suited for heat sealing

## MegaBlock® 0.5/1.2/2.2 ml, PP

- Solvent-resistant, even to DMSO
- Autoclavable
- Square cavities with round bottom wells
- For backup storage of pharmaceutical substances
- For DNA isolation, enzyme and tissue culture applications

## MegaBlock® 1.2 ml, PS transparent

- Transparent for viewing
- Ideal for long-term storage of blood samples
- Pierceable mat cover

## MegaBlock®

	Order No.	Description	Volume/Material	Packaging
<b>NEW</b>	82.1969.002	round, raised wells	0.5 ml PP	56 pcs./case
	82.1970.002	round, raised wells	1.2 ml PS transparent	32 pcs./case
	82.1971.002	round, raised wells	1.2 ml PP	32 pcs./case
	82.1972.002	square wells	2.2 ml PP	32 pcs./case

Sealing foils/mats available

SARSTEDT AG & Co.  
P.O. Box 12 20  
D-51582 Nümbrecht  
Phone +49 2293 305-0  
Fax +49 2293 305-122  
info@sarstedt.com  
www.sarstedt.com