

V-Monovette® Urine

The vacuum system for enclosed urine transfer

NEW!



Comprehensive product range

The V-Monovette® Urine is suited for a variety of applications and available in the size options 100 x 15, 75 x 13 mm and 100 x 13 mm with rounded base and 100 x 15 mm with conical base.

For urine collection, the system includes a 100 ml urine cup and a 3 l urine container with integrated sampling port for enclosed urine transfer to the V-Monovette®. A transfer device for urine sampling from open containers is also available.

Optimal hygiene and convenient handling

An enclosed system significantly improves hygiene and is user friendly. With this system, it is no longer necessary to open a urine cup or container to fill the tube so that contamination with urine is practically excluded.



Ordering Information

| | Order No. | Description | Length/ø in mm | Packaging pcs. inner case/box |
|--|-------------|--|-------------------|----------------------------------|
| Urinalysis | 11.2252.001 | V-Monovette® Urine 4 ml* | 75 x 13 | 50/500 |
| | 11.2352.001 | V-Monovette® Urine 6 ml*, rounded base | 100 x 13 | 50/500 |
| | 11.2452.001 | V-Monovette® Urine 10 ml*, rounded base | 100 x 15 | 50/500 |
| | 11.2552.001 | V-Monovette® Urine 9.5 ml*, conical base | 100 x 15 | 50/500 |
| Urine Culture | 11.2253.001 | V-Monovette® Urine 4 ml with stabilizer*, rounded base | 75 x 13 | 50/500 |
| | 11.2453.001 | V-Monovette® Urine 10 ml with stabilizer*, rounded base | 100 x 15 | 50/500 |
| Accessory | 11.1240 | Transfer device | | 50/500 |
| | 75.562.400 | 100 ml Cup with integrated transfer device, assembled lid, sterile, with warranted sterility due to safety label | | 200/box |
| | 77.575.400 | 3 l Urine container with integrated transfer device | | 30/box |
| Accessory Urine Collection from Catheter | 14.1208.700 | Needle holder | | 100/1.000 |
| | 85.1638.935 | Safety-Multifly® Needle with assembled tubeholder, 21 G x 1¼" tubing 8" | | 100/500 |

* All versions of the V-Monovette® Urine are steril

SARSTEDT AG & Co.
P.O. Box 12 20
D-51582 Nümbrecht
Phone +49 2293 3050
Fax +49 2293 305-2470
info@sarstedt.com
www.sarstedt.com